

## **FACS Analysis of Propidium Iodide Stained Cells to Assess DNA content**

- 1) Spin down  $1 \times 10^7$  cells (about 100 $\mu$ l of an overnight culture) and remove supernatant.
- 2) Add 1 ml of ice cold 70% Ethanol and vortex (cells will appear clumped, this is normal), leave for 10 mins (this stage fixes the cells by dehydration).
- 3) At this stage these cells can be stored at 4<sup>o</sup>C indefinitely if required (good for time courses)
- 4) Spin down cells, remove ethanol solution and re-hydrate by adding 1 ml of 50mM Na Citrate.
- 5) Spin down again, discard supernatant and resuspend pellet in 0.5ml of 50mM Na Citrate containing 1mg/ml RNAase A (as PI binds to all nucleic acids, this step removes the background by digesting RNA in the cell) and incubate for 2h at 37<sup>o</sup>C.
- 6) Add PI to a final concentration of 6 $\mu$ g/ml.
- 7) Load entire contents of the tube into a FACS tube and record fluorescence using the FL2 (red) detector.